


MINI LAPAROSCOPE



Advanced fiber optic technology and improved instrumentation have created a revolution in the field of Minimally Invasive Surgery. This latest development "Mini- or Micro-laparoscopy", uses a mini-laparoscope of only 2 mm in diameter, allowing for much smaller incisions through a Veress needle. In addition to reducing patient trauma, other benefits include reduction of medical costs and greater surgical precision.

The mini-laparoscope can be used for diagnostic laparoscopy, fertility evaluation, the diagnosis of pelvic pain as well as immediate diagnosis in emergency departments.



Features and Benefits

- 30,000 pixel display for greater resolution.
- Semi-flexible shaft improves life-time and increases reliability.
- Small diameter reduces incision size and patient trauma.



CONTEC MEDICAL LTD.
11, Haroshet St.,
P.O. Box 1400
Ramat Hasharon
47113, Israel
Tel: +972-3-549-1153
Fax: +972-3-549-1095
e-mail:
info@contecmedical.com
www.contecmedical.com



MINI LAPAROSCOPE



Technical Specifications

Diameter	1.98 mm (1.7 mm optional)
Working Length	24 cm (9.6 in)
Direction of View	0°
Field of View	80°
Light Cable Connector	Contec/Olympus/Storz/Wolf/ ACMI (all included)

Ordering Information

OP 2002	Semi-Flexible 2 mm Mini Laparoscope, 0°
OP 1603	Mini Cannula for 2 mm Laparoscope, Single-Use, Disposable, (Box of 10)
OP 1602	Disposable Veress Needle, 150 mm (Box of 10)
IS 0015	Reusable Veress Needle, 150 mm (Box of 1)

Hand Instruments



2 mm Semi-Flexible Biopsy Forceps, Serrated Jaws

LI 0200



2 mm Semi-Flexible Biopsy Forceps, Scoop

LI 0201



2 mm Semi-Flexible Grasping Forceps

LI 0202



2 mm Semi-Flexible Scissors

LI 0203

Specifications are subject to change without notice.

Contec Medical designs, manufactures, markets and supports a complete line of high-performance endoscopy products for Minimally Invasive Diagnostic and Surgical procedures, suitable for laparoscopy, arthroscopy, hysteroscopy, thoracoscopy and ENT applications.



CONTEC MEDICAL LTD.

11, Haroshet St.,

P.O. Box 1400

Ramat Hasharon

47113, Israel

Tel: +972-3-549-1153

Fax: +972-3-549-1095

e-mail:

info@contecmedical.com

www.contecmedical.com